

from cylinders

3000 psi max.

Used for connections between insufflators and gas cylinders with post-type valves. Both of these items meet strict Compressed Gas Association (CGA) guidelines.



Pin-indexed Yoke Assembly
CGA 940 x 1/4" SAE Male Flare

part # 9406



36" Flexible Stainless Steel Pigtail
1/4" SAE Female Flare x 1/4" SAE Female Flare

part # FP-123

from outlets

200 psi max.

Used for connections between house gas supply (*i.e.* wall outlets) or regulators and insufflators.



36" Hose Assembly, Conductive, Grey
DISS 1080-A Female x 1/4" SAE Female Flare

This hose assembly is equipped with a one-way check valve to ensure proper use.

part # IN-1080



*Need a different CO2 hose assembly
like DISS female by DISS female?*

Contact Bay Corporation for assistance!!!

 **BAY CORPORATION**[®]
Your Medical Gas Fittings Specialist[™]
867 CANTERBURY ROAD, WESTLAKE, OH 44145-1486 USA

PH: 440.835.2212
FAX: 440.835.1377
888.835.3800
www.baycorporation.com
sales@baycorp.com



INSUFFLATOR Connections