



AY CORPORATION[®]
Your Medical Gas Fittings Specialist[™]

Medical Gas Hose Assemblies



We produce over 100,000 medical gas hose assemblies each year using our full line of medical gas fittings and connections for Oxygen, Medical Air, Nitrous Oxide, Suction (Vacuum), Nitrogen, WAGD (AGSS), Instrument Air, Carbon Dioxide, and Heliox gas services.

Our hose assemblies, which are DEHP free, meet or exceed the specifications designated by the ISO 5359, EN-739 and CSA-Z5359 standards.

We are a trusted partner supplier to most all of the leading global Anesthesia machine and Ventilator original equipment manufacturers, including BD Carefusion, GE Healthcare, Draeger Medical, Hamilton Medical, Medtronic, Philips Healthcare, Smiths Medical and Vyair Medical. In addition, our hose assemblies are used by boom arm, dental sedation, headwall, insufflation and scuba regulator device manufacturers around the world.

Our medical gas hose assemblies are produced utilizing our proprietary medical gas hose, which has been designed, engineered and manufactured to provide a consistent reliable and safe implement for transporting medical gases to its specific function.

We offer a 90-degree swivel option, which was the first one of its kind on the market, for both the inlet and outlet connections.



867 Canterbury Road, Westlake, OH 44145 USA

440-835-2212 sales@baycorp.com